



EVALUATION OF DENTAL RADIOLOGICAL EQUIPMENT
RADIOLOGICAL HEALTH PROGRAM (907) 334-2107
doh.radiation.control@alaska.gov

FORM MUST BE COMPLETED AFTER ANY EQUIPMENT SERVICE OR ASSESSMENT WHICH EFFECTS THE TUBE OR IMAGE QUALITY OF THE DEVICE. Note: This form is intended to document initial installs, annual service, or repairs to a device which involves the tube head. Minor repairs to devices, such as switch replacements, are not to be listed on this form.

PART I. DEVICE STATUS:

INITIAL INSTALL ANNUAL SERVICE REPAIR OTHER _____

PART II. FACILITY INFORMATION:

REGISTRATION NUMBER: _____ Primary dentist license # _____

NAME OF BUSINESS: _____

OR Doing-Business-As: _____

MAILING ADDRESS: _____

CITY/STATE/ZIP: _____

TELEPHONE NUMBER: _____

POINT OF CONTACT: _____

FAX/E-MAIL ADDRESS: _____

PART II. ALL DEVICE(S) MUST BE EVALUATED IN ACCORDANCE TO REGULATIONS (PART F) & MANUFACTURER

RECOMMENDATIONS: The following minimum conditions must be checked for the appropriate device used at a facility (refer to the regulations and/or request a DOH Equipment Inspection Guide). All device parameters must be recorded on a separate document showing each measurement for each device evaluated. Refer to Part III below.

- | | | | |
|--|---------------------------|---|--|
| ✓ Device is FDA registered/permitted | ✓ Provider qualified | ✓ Quality assurance measurements | ✓ Warning label |
| ✓ Radiation leakage from source tube | ✓ Air kerma emitted | ✓ Technique factor | ✓ Beam quality (HVL) |
| ✓ Aluminum equivalent measurement | ✓ Fluoroscopic filtration | ✓ Battery charge indicator | ✓ Modified components |
| ✓ Mechanical support of tube head | ✓ Multiple tubes | ✓ Source to skin distance (SID <2%) | ✓ Locks |
| ✓ Radiation exposure control | ✓ Exposure initiation | ✓ Exposure position | ✓ Dark room |
| ✓ Backscatter shields on handheld | ✓ Beam on indicators | ✓ Exposure reproducibility (<5%) | ✓ Timer- Rotation check |
| ✓ Kilovolt peak accuracy (<10%/<20%) | ✓ Beam alignment (CT) | ✓ Device accreditation (CT) | ✓ Shutter check (CT) |
| ✓ Termination of exposure (timer/switch) | ✓ Tomographic plane (CT) | ✓ Other checks as regulations or device manufacture require | ✓ Device pass /fail rating coefficient of variance |
| | ✓ Intraoral used kVp >51 | | |

PART III. RADIATION PRODUCING EQUIPMENT: Devices requiring installation, evaluation, or preventative maintenance must include the following device identification information and all applicable measurements described in Part II above: Manufacture, Model, Serial Number, and Device Type (such as CT, CBCT, PANO, intra-oral, etc.) **DEVICE PARAMETERS MUST BE RECORDED ON A SEPERATE DOCUMENT SHOWING EACH MEASUREMENT FOR EACH DEVICE EVALUATED. EACH DEVICE MUST STATE A PASS OR FAIL RATING AND ANY CORRECTIVE ACTIONS, IF FAILED. ATTACH ALL DOCUMENTS RECORDING MEASUREMENTS TO THIS FORM.**

PART IV. CERTIFICATION OF QUALIFIED SERVICE PROVIDER/VENDOR/OR INSPECTOR: THIS IS TO CERTIFY THAT, I, THE DELEGATED AUTHORITY, TO THE BEST OF MY KNOWLEDGE AND BELIEF, ALL INFORMATION CONTAINED HEREIN, INCLUDING ANY SUPPLEMENTS ATTACHED HERETO, IS CORRECT AND MEETS THE REQUIREMENTS OF THE SUGGESTED STATE REGULATIONS-PART F- PUBLISHED BY THE CONFERENCE OF RADIATION CONTROL PROGRAM DIRECTORS (CRCPD).

PRINT NAME: _____ SIGNATURE: _____

COMPANY NAME _____ Date: _____

COMPANY EMAIL ADDRESS & PHONE # _____